Hiden Analytical Ltd. 420 Europa Boulevard Warrington WA5 7UN England T +44 [0] 1925 445 225
F +44 [0] 1925 416 518
E info@hiden.co.uk
W www.HidenAnalytical.com

Gas Analysis Performance Data Sheet PDS-300012

## **Results for a Standard QGA Pre and Post Bakeout**

m/z	Species of Interest	Pre-bake ppm	Post-bake ppm
Carrier	<sup>40</sup> Ar	100%	100%
2	H <sub>2</sub>	1147.9	1508.1
3	HD/H₃	0.2	0.3
4	He	<0.1	<0.1
12	С	1.2	0.7
13	<sup>13</sup> C	8.3	8.6
14	N	1.3	0.9
15	CH <sub>3</sub>	2.1	1.9
16	CH <sub>4</sub> /O	8.3	3.9
17	NH <sub>3</sub> /OH	65.5	11.8
18	H <sub>2</sub> O/ <sup>36</sup> Ar <sup>++</sup>	1024.3	750.6
18	H <sub>2</sub> O	284.6	51.2
19	H <sub>3</sub> O/F	159.7	154.6
20	Ne/Ar <sup>++</sup>	186845	182088
28	N <sub>2</sub> /CO	26.7	21.7
30	NO	1.6	0.4
32	$O_2$	3.5	1.1
34	PH <sub>3</sub> /H <sub>2</sub> S	<0.1	<0.1
35	<sup>35</sup> Cl	6.6	10.1
36	<sup>36</sup> Ar	3525.9	3366.5
37	<sup>37</sup> Cl	2.2	2.9
38	<sup>38</sup> Ar	652.5	610.8
44	CO <sub>2</sub>	16.9	7.7
45	<sup>13</sup> CO <sub>2</sub>	0.3	0.1
57	Hydrocarbons	<0.1	<0.1
64	<sup>32</sup> SO <sub>2</sub>	0.1	<0.1
69	CF <sub>4</sub>	<0.1	<0.1
78	C <sub>6</sub> H <sub>6</sub>	<0.1	<0.1
84	Kr	<0.1	<0.1



The table of results shows data collected post-bakeout for a standard QGA system. This system consisted of:

- HAL 201 gauge with an open ion source
- Standard capillary
- EXT75H turbo
- nXDS6i scroll pump

The system was run on BIP+ Argon for a period of 72 hours, all impurities in the carrier gas are <0.01 ppm. The system manifold then underwent an overnight (16 hours) bake-out and the measurements repeated over another 72-hour period. During testing, total pressure of the system was 3.5E-6 mbar.